

QuikChange Site-Directed Mutagenesis

Please refer to Maniatis for theoretical background. This rapid four-step procedure generates mutants with greater than 80% efficiency. The protocol is simple and uses miniprep plasmid DNA. The mutagenesis protocol is used to make point mutations, switch amino acids, and delete or insert single or multiple amino acids. The QuikChange site-directed mutagenesis method is performed using a *proof-reading* DNA polymerase cyclase. The basic procedure utilizes a supercoiled double-stranded DNA (dsDNA) vector with an insert of interest and two synthetic oligonucleotide primers containing the desired mutation. The oligonucleotide primers, each complementary to opposite strands of the vector, are extended during temperature cycling by DNA polymerase. Incorporation of the oligonucleotide primers generates a mutated plasmid containing staggered nicks. Following temperature cycling, the product is treated with *Dpn* I. The *Dpn* I endonuclease (target sequence: 5'-Gm6ATC-3') is specific for methylated and hemimethylated DNA and is used to digest the parental DNA template and to select for mutation-containing synthesized DNA. DNA isolated from almost all *E. coli* strains is *dam* methylated and therefore susceptible to *Dpn* I digestion. The nicked vector DNA containing the desired mutations is then transformed into XL1-Blue competent cells.


Flow chart in vitro mutagenesis

The protocol depends on PCR, therefore the reaction itself is bound to create new mutations. Therefore, subclone the fragment in which mutations should be introduced into Bluescript or related general cloning vector. For this purpose, use unique sites.

Perform mutagenesis as described and analyse colonies by restriction analysis. Choose 1 or 2 clones for sequencing the insert. If sequence is correct, cut out fragment with same enzymes as before and re-insert into original vector.

Prepare minipreps, check by restriction (same sites as above plus 1 site in fragment, ideally one that was introduced with mutation). If correct, go ahead with experiment.

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