

PCR

Primer

The most important thing to establish a new PCR is to choose your primers. There are programs in the net like Primer 3 (http://frodo.wi.mit.edu/cgi-bin/primer3/primer3_www.cgi) that do this for you however, primers should be between 15 and 25 bp long with an optimum at 20 bp. They should anneal between 55 and 65°C with an optimum at 60°C, and they should have an GC percentage of about 50%.

Primers are delivered lyophilized and with tare solubilized in Sigma water as 100µM Stock. Aliquot the primers from the stock as to reduce repeated thawing of the stock. The primers are used at a 25mM concentration in the mix, so always freshly dilute the primers.

PCR Setup

All ingredients are to be equilibrated to 37°C, vortexed, spun down and directly put on ice. The entire PCR is pipetted on ice!!!! Especially the Taq polymerase has to be kept ice cold at all times and therefore is pipetted last!!!

Ansatz:

SAMPLE NUMBER COMPONENT	1X ANSATZ [ML]
Leber cDNA	
10 x Buffer	2,5
50 mM MgCl₂	0,75
5 mM dNTPs	1
Primer F [25µM]	1
Primer R [25µM]	1
DMSO	
Taq Pol	0,2
H₂O up to 25µl	17,55

PCR condition depend on size of the fragment, primers and DNA used.

Taq polymerase needs 30'' to synthesize 500bp. So set your elongation time appropriately. The best annealing temperature for the primers has to be tested over a temperature gradient. For genomic DNA always use 35 cycles, for Plasmid DNA 15 cycles are sufficient.

PCR condition - example:

Denature	94°C 2'
Denature	94°C 30''
Anneal	65°C 30''
Elongate	72°C 90''
Elongate	72°C 5'

PCR conditions can vary, if the initial temperature gradient does not work you can additionally titrate various compounds: 1-5µl DNA; 1- 2.5mM MgCl₂, 0-10% DMSO. Have a look at Henry Sun's paper (*Tung-Tien Sun, Excessive trust in authorities and its influence on*

