

Addition of 3' A- overhangs to PCR products after amplification

Introduction

Direct cloning of DNA amplified by proofreading polymerases into pET SUMO is often difficult because proofreading polymerases remove the 3' A-overhangs necessary for TA Cloning®. Invitrogen has developed a simple method to clone these blunt-ended fragments.

Before Starting

You will need the following items:

- *Taq* polymerase
 - A heat block equilibrated to 72°C
 - Phenol-chloroform (optional)
 - 3 M sodium acetate (optional)
 - 100% ethanol (optional)
 - 80% ethanol (optional)
 - TE buffer (optional)
-

Procedure

This is just one method for adding 3' adenines. Other protocols may be suitable.

1. After amplification with a proofreading polymerase, place vials on ice and add 0.7-1 unit of *Taq* polymerase per tube. Mix well. It is not necessary to change the buffer. A sufficient number of PCR products will retain the 3' A-overhangs.
2. Incubate at 72°C for 8-10 minutes (do not cycle).
3. Place on ice and use immediately in the ligation reaction.

Note: If you plan to store your sample overnight before proceeding with the ligation reaction, extract your sample with an equal volume of phenol-chloroform to remove the polymerases. Ethanol-precipitate the DNA and resuspend in TE buffer using the starting volume of the PCR.



Note

You may also gel-purify your PCR product after amplification with a proofreading polymerase. After purification, add *Taq* polymerase buffer, dATP, and 0.5 unit of *Taq* polymerase. Incubate the reaction for 10-15 minutes at 72°C and use in the ligation reaction.
