## Isolation of total RNA

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### **Introduction:**

**Yields to expect: 1-3 µg of total RNA/mg of tissue. Keep in mind that mRNA represents 1-3 % of total RNA. Methods 1 and 2 are equivalent.**

**For standard preparations of RNA, from which you expect <1.000 µg of total RNA, use method 1. If larger amounts of total RNA are required, you must use protocol 3 (one mouse liver for example gives 4-5 mg of total RNA.**

**For most applications (cDNA synthesis, Northern blotting, hybrid selection), total RNA is fine. In cases where large amounts of mRNA are needed, enrichment by oligo dT selection (protocol 4) is needed.**

**Whatever method you use, RNA is always stored as EtOH precipitate. This means that the final RNA, resuspended in ultrapure Sigma water, is reppt. with 0.1 vol. of 3 M NaAc, 2 vol. EtOH and stored at -20°C. If RNA is needed, vortex this suspension and remove appropriate amount in fresh tube. Then spin as usual, wash with 70 % EtOH, air-dry and resuspend in Sigma water for further use.**

### **1. TRIZOL<sup>®</sup> Reagent**

The protocol should not be used for large numbers of cells (more than 50 Mio.) or more than 1 g of tissue. In this case, refer to “stepwise RNA preparation protocol!

#### **1.1 Homogenization:**

##### **a. Tissue:**

Homogenize tissue samples in at least **1 mL TRIZol Reagent per 50-100 mg** of tissue using a **Polytron** for 2 min at full speed (cool intermittently). *Polytron: Use safety goggles! Before and*

*after use: run for 1 min in 10 N NaOH, followed by 2-3 x washes with sterile water. Finally, rinse with EtOH.* The sample volume should not exceed 10% of the volume of TRIzol reagent used for homogenization. It is always better to use more Trizol in order to improve yield and purity of RNA. **Always add Vanadyl-Ribonucleotide** (200 mM stock, Invitrogen) at 5 mM final concentration as additional protection against RNase.

#### **b. Cells grown in monolayer:**

Wash cells with sterile PBS. Lyse cells directly in a culture dish by adding 1 mL TRIzol reagent per 10 cm<sup>2</sup>. (*The amount of TRIzol reagent added is based on the area of the culture dish, at least →1mL / 10 cm<sup>2</sup>. An insufficient amount of TRIzol reagent results in contamination with DNA and low yield.*)

Pass the cell lysate several times through a pipette and use polytron instantly. Never trypsinize cells for RNA prep.! If you need to collect cells, scrape them.

#### **c. Cells grown in suspension:**

Pellet cells by centrifugation.

Wash cells with sterile PBS before adding TRIzol. Use 1 mL of TRIzol Reagent per 5-10 x 10<sup>6</sup> of animal or plant cells; 1 x 10<sup>7</sup> bacterial cells. Pass the cell lysate several times through a pipette and use polytron instantly.

#### **Optional:**

An additional isolation step may be required for samples with high content of protein. Remove insoluble material by centrifuging at 12.000 x g for 10 min. at 4°C. The supernate contains the RNA. Transfer the cleared homogenate to a fresh tube. Continue with phase separation step.

### **1.2 Phase separation:**

- **Incubate** the homogenized samples for **5 min at 15-30°C** to permit complete dissociation of nucleoprotein complexes.
- **Add 0,2 mL of chloroform** (per 1 mL TRIzol). Close lid. Shake tubes vigorously by hand for approx. 15 sec. and incubate at 15-30°C for 2-3 min.
- **Centrifuge** samples at 12.000 x g for 10 min. at 4°C → Phase separation appears.
  - lower red Phase = phenol-chloroform (organic)
  - interphase

colorless upper aqueous phase = RNA

- Carefully transfer all aqueous phase into fresh tube.

### 1.3 RNA precipitation:

- **Transfer** upper phase to a fresh tube (save organic phase if isolation of DNA / protein is desired).
- Mix with **Isopropanol: 0,5 mL** isopropanol (per 1 mL TRIzol). Incubate at RT for 10 min → centrifuge at 12.000 x g for 10 min at 4°C.
- RNA precipitate forms a gel-like pellet on side and bottom.

### 1.4 RNA Wash: (~EtOH preequilibrated at 4°C)

- Remove supernatant
- Wash RNA pellet 1x with 75% EtOH, (at least 1 mL per 1mL TRIzol)
- Mix by vortexing
- centrifuge at no more than 7.500 x g for 2 min. 4°C.
- Air-dry the pellet

### 1.5 Redissolving the RNA:

Dissolve RNA in appropriate vol. (100- 1.000 µl) of RNase free water („RNA grade“) by passing a few times through pipette and incubating for 5-10 min. at 55-60°C. Reextract at least 1 x with Ph/Chl: Add RNA grade water to 500-1.000 µl and 1 µl of glycogen. Add 0.8 vol. of Ph/Chl, vortex 10 sec. And spin 5 min at 12.000 x g. Transfer supernatant in fresh tube. If there's still interphase or dark colour (from Vanadyl-ribonucleotide), repeat this step.

Then, add 1 vol. of Chl., vortex and spin as above. Transfer supernatant in fresh tube. Add 0.1 vol. of 3 M NaAc and 2 vol of EtOH. Mix and spin 10 min at 12.000x g. Discard supernatant and wash RNA pellet with 0.5 ml of 70 % EtOH. Vortex, spin briefly and remove spnt completely. Air-dry. Resuspend in 100-500 µl of RNA grade water.

### 1.6 Photometry and storage:

Determine amount by **photometry** (use microcuvettes, Photometry in TE ( $1 \text{ OD}_{260\text{nm}} = 40 \text{ }\mu\text{g/ml}$  RNA; OD 260/280 should be more than 2) and quality by electrophoresis (Load 1, 2 and 5  $\mu\text{l}$ ).

**For storage**, add 0.1 vol of 3 M NaAc pH 5.5, 2 vol. EtOH, mix and keep at  $-20^\circ\text{C}$ .

**Never store RNA in water!**

To remove a given amount of RNA for experiments: vortex briefly, remove desired amount of suspension and transfer in Eppendorf tube. If working with less than 10  $\mu\text{g}$ , add 1  $\mu\text{l}$  glycogen. Spin 10 min at 14.000 rpm at RT, remove spnt. Wash 1 x with 70 % EtOH, spin briefly, discard and air-dry RNA pellet.

**Remarks:**

1. Isolation of RNA from small quantities of tissue  $\rightarrow < 1$  to 10 mg or cells  $< 10^2$  to  $10^4$  add 800  $\mu\text{l}$  of TRIzol Reagent to the tissue or cells !

Add 20  $\mu\text{g}$  glycogen with final concentration in TRIzol of 250  $\mu\text{g/mL}$ . The glycogen remains in the upper aqueous phase.

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