

# Protocol

Theme: Double transfection of Keratinocytes with Lipofectamine™ 2000 Invitrogen

Equipent: Lipofectamin™ 2000 Invitrogen, DNA/plasmid of interest, medium (full medium, reduced medium, minimal medium), 24-well plate

Procedure: for 24-well format, all amounts and volumes are given on a per well basis, composition of medium is given on the next page

1. 1 day before transfection, plate  $1,2 \times 10^5$  cells/well in 500  $\mu$ l **reduced medium** and incubate at 32°C and 5 % CO<sub>2</sub>.
  2. Next day discard the old medium, wash twice with **minimal medium** and add 500  $\mu$ l **minimal medium** to each well.
  3. **For each transfection sample**, prepare complexes as follows:
    - a. Dilute **1,6  $\mu$ g precipitated DNA** in 50  $\mu$ l **minimal medium** and mix gently.
    - b. Mix Lipofectamine™ 2000 gently before use, then dilute **2  $\mu$ l Lipofectamin** in 50  $\mu$ l of **minimal medium**. Incubate for 5 minutes at room temperature.
    - c. After the 5 minute incubation, combine the diluted DNA with diluted Lipofectamine™ 2000 (total volume = 100  $\mu$ l). Mix gently and incubate for 20 minutes at room temperature (solution may appear cloudy).
- Note:** Complexes are stable for 6 hours at room temperature.
4. Add the 100  $\mu$ l of complexes to each well containing cells and medium. Mix gently by rocking the plate back and forth.
  5. Incubate cells at 32°C and 5 % CO<sub>2</sub> for 6 hours.
  6. Change medium after 6 hours to **reduced medium**.
  7. Incubate cells at 32°C and 5 % CO<sub>2</sub> over night.
  8. Next day retransfect the cells as described above (**Step 2-5**).
  9. After the second transfection change medium to **full medium** and incubate the cells at 32°C and 5 % CO<sub>2</sub> prior to testing for transgene expression (24-72 hours).

With this procedure you can expect a transfection efficacy of ~ 40%.

## Scaling Up or Down Transfections

Culture vessel	Surf. area per well <sup>1</sup>	Shared reagents		DNA transfection	
		Vol. of plating medium	Vol. of dilution medium <sup>2</sup>	DNA	Lipofectamin 2000
96-well	0.3 cm <sup>2</sup>	100 µl	2 x 25 µl	0.4 µg	0.5 µl
24-well	2 cm <sup>2</sup>	500 µl	2 x 50 µl	1.6 µg	2.0 µl
12-well	4 cm <sup>2</sup>	1 ml	2x 100 µl	3.2 µg	4.0 µl
6-well	10 cm <sup>2</sup>	2 ml	2 x 250 µl	8 µg	10 µl
60-mm	20 cm <sup>2</sup>	5 ml	2 x 0.5 ml	16 µg	20 µl
10-cm	60 cm <sup>2</sup>	15 ml	2 x 1.5 ml	48 µg	60 µl

<sup>1</sup> Surface areas may vary depending on the manufacturer.

<sup>2</sup> Volumes of dilution medium in Step 3a & 3b.

**Full Medium:**

460 ml	DMEM low Calcium (50µM)
50 ml	FCS (Chelex treated)
5 ml	Glutamax (100x)
5 ml	Pyruvat (100x)
2,5 ml	Pen/Strep
2 ml	Adenin (250x) 45 mM
500 µl	EGF (1000x) 10µg/ml
500 µl	Insulin (1000x) 5 mg/ml
250 µl	Hydroxcortison (2000x)
5 µl	Choleratoxin (10 <sup>-5</sup> M)

**Reduced Medium:**  
(without Pen/Strep)

460 ml	DMEM low Calcium (50µM)
50 ml	FCS (Chelex treated)
5 ml	Glutamax (100x)
5 ml	Pyruvat (100x)
2 ml	Adenin (250x) 45 mM
500 µl	EGF (1000x) 10µg/ml
500 µl	Insulin (1000x) 5 mg/ml
250 µl	Hydroxcortison (2000x)
5 µl	Choleratoxin (10 <sup>-5</sup> M)

**Minimal Medium:**  
(without Pen/Strep and FCS)

460 ml	DMEM low Calcium (50µM)
5 ml	Glutamax (100x)
5 ml	Pyruvat (100x)
2 ml	Adenin (250x) 45 mM
500 µl	EGF (1000x) 10µg/ml
500 µl	Insulin (1000x) 5 mg/ml
250 µl	Hydroxcortison (2000x)
5 µl	Choleratoxin (10 <sup>-5</sup> M)